


## GLP-1R agonist 8

## Chemical Properties

CAS No. :	2736446-82-3
Formula:	C33H32N4O5
Molecular Weight:	564.63
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	GLP-1R agonist 8 is a potent GLP-1R agonist with an EC50 of less than 2 nM.
Targets(IC50)	Others,Glucagon Receptor

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7711 mL	8.8554 mL	17.7107 mL
5 mM	0.3542 mL	1.7711 mL	3.5421 mL
10 mM	0.1771 mL	0.8855 mL	1.7711 mL
50 mM	0.0354 mL	0.1771 mL	0.3542 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481