

Compound N016-0043

## Chemical Properties


CAS No. :

Formula: C<sub>60</sub>H<sub>95</sub>N<sub>3</sub>O<sub>19</sub>S<sub>3</sub>

Molecular Weight: 1258.6

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Compound N016-0043 is a useful organic compound for research related to life sciences and the catalog number is T131340.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7945 mL	3.9727 mL	7.9453 mL
5 mM	0.1589 mL	0.7945 mL	1.5891 mL
10 mM	0.0795 mL	0.3973 mL	0.7945 mL
50 mM	0.0159 mL	0.0795 mL	0.1589 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481