

Fmoc-Ala-Ala-Asn-PABC-PNP

Chemical Properties

CAS No. : 1834516-06-1


Formula: C₃₉H₃₈N₆O₁₁

Molecular Weight: 766.75

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Fmoc-Ala-Ala-Asn-PABC-PNP is a cleavable peptide ADC linker [1].
Targets(IC50)	ADC Linker,PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3042 mL	6.521 mL	13.0421 mL
5 mM	0.2608 mL	1.3042 mL	2.6084 mL
10 mM	0.1304 mL	0.6521 mL	1.3042 mL
50 mM	0.0261 mL	0.1304 mL	0.2608 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Shaosong Chu. Peptide-drug conjugates. US20150343083A1.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481