


Tivantinib

Chemical Properties

CAS No. :	905854-02-6
Formula:	C ₂₃ H ₁₉ N ₃ O ₂
Molecular Weight:	369.42
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Tivantinib (ARQ 197) is an orally bioavailable small molecule inhibitor of c-Met with potential antineoplastic activity.
Targets(IC50)	Apoptosis,c-Met/HGFR
In vitro	Treatment with ARQ-197 across three heterotrophic transplantation models demonstrated a reduction in tumor growth: 66% in the HT29 model, 45% in the MKN-45 model, and 79% in the MDA-MB-231 model. Oral administration of 200 mg/kg ARQ-197 in these HT29, MKN-45, and MDA-MB-231 xenograft models did not result in significant weight changes. Pharmacologically, ARQ-197 notably inhibited the phosphorylation of c-Met in the human colon xenograft tumor HT29. A single oral dose of 200 mg/kg ARQ-197 led to a strong decrease in c-Met auto-phosphorylation after 24 hours. In summary, ARQ-197 inhibits the growth of c-Met-dependent human tumor xenografts.
In vivo	ARQ-197 exhibits antitumor activity by inhibiting the proliferation of A549, DBTRG, and NCI-H441 cells, with IC50 values of 0.38, 0.45, and 0.29 μ M respectively. Treatment with ARQ-197 reduces phosphorylation in the MAPK signaling cascade, thereby preventing invasion and migration. Moreover, ectopic expression of c-Met in the NCI-H661 cell line, which lacks endogenous c-Met expression, leads to an invasive phenotype that is inhibited by ARQ-197. This compound also inhibits human recombinant c-Met with a constant K_i value of 355 nM and blocks downstream c-Met signaling pathways by inhibiting c-Met phosphorylation. Despite increased concentrations of ARQ-197 not significantly affecting the K_m of ATP, exposure to 0.5 μ M ARQ-197 results in a threefold decrease in the V_{max} of c-Met. The decrease in V_{max} without affecting ATP's K_m suggests ARQ-197 inhibits c-Met through a non-ATP competitive mechanism, indicating high kinase selectivity. ARQ-197 suppresses both constitutive and ligand-mediated c-Met autophosphorylation, thereby inhibiting c-Met activity and downstream effectors. In human cancer cells expressing c-Met, including HT29, MKN-45, and MDA-MB-231, ARQ-197 induces an increase in caspase-dependent apoptosis.
Kinase Assay	c-Met SDS-PAGE in vitro kinase assay: Recombinant c-Met protein (100 ng) is preincubated with increasing concentrations of ARQ-197 for 30 minutes at room temperature. Following preincubation, 100 μ M of poly-Glu-Tyr substrate and various concentrations of ATP containing 5 μ Ci of [γ - ³² P]ATP are added to the reaction mixture. The reaction is incubated for 5 minutes at room temperature and then stopped by the addition of 5 μ L of SDS-polyacrylamide gel, reducing sample buffer. The samples are

Kinase Assay	then loaded onto a 7.5% acrylamide gel and SDS-PAGE is performed. The phosphorylated poly-Glu-Tyr substrates are ultimately visualized by autoradiography. c-Met activity is quantified by densitometry.
Cell Research	HT29, MKN-45, and MDA-MB-231 cells are seeded in black 96-well plates at 5×10^3 cells per well overnight in a medium with 10% FBS. The next day, cells are treated with increasing concentrations of ARQ-197 (0.03-10 μ M) for 24, 32, and 48 hours at 37 °C. After ARQ-197 treatment, the drug-containing medium is removed and cells are incubated for at least 10 minutes in a labeling solution (10 mM HEPES, 140 mM NaCl, and 6 mM CaCl ₂) containing 2 μ g/mL Hoescht 33342 (blue channel), 500-times diluted Annexin V-FITC (green channel), and 1 μ g/mL propidium iodide (red channel). High-content image acquisition and analysis are carried out. The program is set to take four images per well. The exposure time is set at 16.7 ms/10% gain, 500 ms/35% gain, and 300 ms/30% gain for the 4,6-diamidino-2-phenylindole, FITC, and rhodamine channels, respectively. Images are processed and the numbers of positive cells for each channel and each condition are determined. In addition, HT29 cells are treated with increasing concentrations of ARQ-197 for 32 hours in the absence or the presence of 25, 50, and 100 μ M ZVAD-FMK (irreversible general caspase inhibitor), and the same procedures are undertaken. All experiments are done in triplicate. To determine whether the apoptotic effect is due to c-Met inhibition, the effect of ARQ-197 when glyceraldehyde-3-phosphate dehydrogenase (GAPDH) and c-Met are knocked down using siRNA is investigated. HT29, MKN-45, and MDA-MB-231 cells are transfected with a nontargeted control siRNA, a gapgh-targeted control siRNA, or a met-targeted siRNA. After 3 days, c-Met, GAPDH, and β -actin expression levels are determined using specific antibodies. To determine if the effect is caspase dependent, HT29, MKN-45, and MDA-MB-231 cells are transfected with a met-targeted siRNA for 2 days and incubated in the absence or the presence of increasing concentrations of ZVAD-FMK for 1 additional day. A nontargeted siRNA and a gapgh-targeted siRNA (siRNA GAPDH) are also transfected in parallel, as controls. Cells are then stained with Annexin V-FITC and propidium iodide, and the percentage of apoptotic cells is determined. (Only for Reference)

Solubility Information

Solubility	DMSO: 125 mg/mL (338.37 mM), Sonication is recommended. Ethanol: < 1 mg/mL (insoluble or slightly soluble), H ₂ O: < 1 mg/mL (insoluble or slightly soluble), (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (5.41 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7069 mL	13.5347 mL	27.0695 mL
5 mM	0.5414 mL	2.7069 mL	5.4139 mL
10 mM	0.2707 mL	1.3535 mL	2.7069 mL
50 mM	0.0541 mL	0.2707 mL	0.5414 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Munshi N, et al. Mol Cancer Ther. 2010, 9(6), 1544-1553.

Huang Y, Guo Y, Zhou Y, et al. Tivantinib alleviates inflammatory diseases by directly targeting NLRP3. iScience. 2023

Comoglio PM, et al. Nat Rev Drug Discov, 2008, 7(6), 504-516.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481