


## Xylamidine tosylate

## Chemical Properties

CAS No. : 6443-40-9  
 Formula: C<sub>26</sub>H<sub>32</sub>N<sub>2</sub>O<sub>5</sub>  
 Molecular Weight: 484.61  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Xylamidine tosylate is a Serotonin inhibitor with reproductive effects
Targets(IC50)	Others,Endogenous Metabolite

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0635 mL	10.3176 mL	20.6351 mL
5 mM	0.4127 mL	2.0635 mL	4.127 mL
10 mM	0.2064 mL	1.0318 mL	2.0635 mL
50 mM	0.0413 mL	0.2064 mL	0.4127 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481