


3,5-Di(prop-1-ynyl)pyridine

Chemical Properties

CAS No. :	1822984-25-7
Formula:	C ₁₁ H ₉ N
Molecular Weight:	155.2
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Heterocyclic compound - pyridine
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.4433 mL	32.2165 mL	64.433 mL
5 mM	1.2887 mL	6.4433 mL	12.8866 mL
10 mM	0.6443 mL	3.2216 mL	6.4433 mL
50 mM	0.1289 mL	0.6443 mL	1.2887 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481