

Compound T131045

Chemical Properties


CAS No. :

Formula: C₃₆H₃₆O₁₈

Molecular Weight: 756.67

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

(3R,4S,5S,6R)-4,5-dihydroxy-2-[[5-hydroxy-2-(4-hydroxyphenyl)-4-oxo-7-[[[(2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]-4H-chromen-3-yl]oxy]-6-(hydroxymethyl)oxan-3-yl (2E)-3-(4-hydroxyphenyl)prop-2-enoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3216 mL	6.6079 mL	13.2158 mL
5 mM	0.2643 mL	1.3216 mL	2.6432 mL
10 mM	0.1322 mL	0.6608 mL	1.3216 mL
50 mM	0.0264 mL	0.1322 mL	0.2643 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481