


2-(Tetrahydro-pyran-4-yl)-ethanol

Chemical Properties

CAS No. :	4677-18-3
Formula:	C7H14O2
Molecular Weight:	130.19
Storage:	Pure form: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-(Tetrahydro-pyran-4-yl)-ethanol, with CAS No. 4677-18-3, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 2-(Tetrahydro-pyran-4-yl)-ethanol provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.6811 mL	38.4054 mL	76.8108 mL
5 mM	1.5362 mL	7.6811 mL	15.3622 mL
10 mM	0.7681 mL	3.8405 mL	7.6811 mL
50 mM	0.1536 mL	0.7681 mL	1.5362 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481