

## Vineridine

## Chemical Properties


CAS No. :

Formula: C<sub>22</sub>H<sub>26</sub>N<sub>2</sub>O<sub>5</sub>

Molecular Weight: 398.459

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Vineridine is a natural product that can be used as a reference standard.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5097 mL	12.5483 mL	25.0966 mL
5 mM	0.5019 mL	2.5097 mL	5.0193 mL
10 mM	0.251 mL	1.2548 mL	2.5097 mL
50 mM	0.0502 mL	0.251 mL	0.5019 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481