


Lipiferolide

Chemical Properties

CAS No. :	41059-80-7
Formula:	C ₁₇ H ₂₂ O ₅
Molecular Weight:	306.35
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Lipiferolide is an inhibitor of farnesyl protein transferase(FPTase), it inhibits the FPTase activity in a dose-dependent manner, and shows cell growth inhibitory activity against several tumor cells, it possesses cytotoxic activity against KB cells. Liriodenine also has antiplasmodial activity.
Targets(IC50)	Parasite,Antifection,Telomerase,Transferase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2642 mL	16.3212 mL	32.6424 mL
5 mM	0.6528 mL	3.2642 mL	6.5285 mL
10 mM	0.3264 mL	1.6321 mL	3.2642 mL
50 mM	0.0653 mL	0.3264 mL	0.6528 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481