


## 4-O-Demethylisokadsurenin D

## Chemical Properties

CAS No. :	89104-59-6
Formula:	C <sub>20</sub> H <sub>22</sub> O <sub>5</sub>
Molecular Weight:	342.391
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	4-O-Demethylisokadsurenin D is a natural product for research related to life sciences. The catalog number is TN3054 and the CAS number is 89104-59-6.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9206 mL	14.6032 mL	29.2065 mL
5 mM	0.5841 mL	2.9206 mL	5.8413 mL
10 mM	0.2921 mL	1.4603 mL	2.9206 mL
50 mM	0.0584 mL	0.2921 mL	0.5841 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481