


Y-27152

Chemical Properties

CAS No. : 127408-30-4
 Formula: C₂₁H₂₂N₂O₄
 Molecular Weight: 366.41
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Y-27152 is a prodrug of the KATP (Kir6) channel opener Y-26763 and is a long-acting K ⁺ channel opener.
Targets(IC50)	Others,Potassium Channel

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7292 mL	13.6459 mL	27.2918 mL
5 mM	0.5458 mL	2.7292 mL	5.4584 mL
10 mM	0.2729 mL	1.3646 mL	2.7292 mL
50 mM	0.0546 mL	0.2729 mL	0.5458 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481