


## N2-(Isopropylphenoxyacetyl)-2'-O-propargylguanosine

## Chemical Properties

CAS No. :	2095417-23-3
Formula:	C <sub>24</sub> H <sub>27</sub> N <sub>5</sub> O <sub>7</sub>
Molecular Weight:	497.5
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Nucleoside Derivatives - 2'-Modified Nucleosides
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0101 mL	10.0503 mL	20.1005 mL
5 mM	0.402 mL	2.0101 mL	4.0201 mL
10 mM	0.201 mL	1.005 mL	2.0101 mL
50 mM	0.0402 mL	0.201 mL	0.402 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481