


Chlortoluron-mecoprop mixt.

Chemical Properties

CAS No. :	37341-10-9
Formula:	C ₂₀ H ₂₄ Cl ₂ N ₂ O ₄
Molecular Weight:	427.32
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Chlortoluron-mecoprop mixt. is a bioactive chemical mixture.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3402 mL	11.7008 mL	23.4017 mL
5 mM	0.468 mL	2.3402 mL	4.6803 mL
10 mM	0.234 mL	1.1701 mL	2.3402 mL
50 mM	0.0468 mL	0.234 mL	0.468 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Tsvetkova OA, Kolosova KY. [Exacerbation of chronic obstructive pulmonary disease: the role of the molecular genetic method for the identification of pathogenic agents]. *Klin Med (Mosk)*. 2016;94(10):770-4. Russian. PubMed PMID: 30299032.

Drugs and Lactation Database (LactMed) [Internet]. Bethesda (MD): National Library of Medicine (US); 2006-. Available from <http://www.ncbi.nlm.nih.gov/books/NBK501887/> PubMed PMID: 30000947.

Drugs and Lactation Database (LactMed) [Internet]. Bethesda (MD): National Library of Medicine (US); 2006-. Available from <http://www.ncbi.nlm.nih.gov/books/NBK501799/> PubMed PMID: 30000858.

Drugs and Lactation Database (LactMed) [Internet]. Bethesda (MD): National Library of Medicine (US); 2006-. Available from <http://www.ncbi.nlm.nih.gov/books/NBK501738/> PubMed PMID: 30000797.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481