


Schidigerasaponin F1

Chemical Properties

CAS No. :	266998-42-9
Formula:	C44H72O18
Molecular Weight:	889.042
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Schidigerasaponin F1 is a useful organic compound for research related to life sciences. The catalog number is T125155 and the CAS number is 266998-42-9.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1248 mL	5.624 mL	11.2481 mL
5 mM	0.225 mL	1.1248 mL	2.2496 mL
10 mM	0.1125 mL	0.5624 mL	1.1248 mL
50 mM	0.0225 mL	0.1125 mL	0.225 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481