


Salicylaldoxime

Chemical Properties

CAS No. :	94-67-7
Formula:	C7H7NO2
Molecular Weight:	137.14
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Salicylaldoxime is a chelating agent sometimes used in the analysis of samples containing transition metal ions, with which it often forms brightly colored coordination complexes.
Targets(IC50)	Others
In vivo	Hhaemolysis of sheep erythrocytes by antibody and either guinea pig or human complement is inhibited by salicylaldoxime (O-hydroxybenzaldoxime.)

Solubility Information

Solubility	DMSO: 55 mg/mL (401.05 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween-80+45% Saline: 5 mg/mL (36.46 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.2918 mL	36.4591 mL	72.9182 mL
5 mM	1.4584 mL	7.2918 mL	14.5836 mL
10 mM	0.7292 mL	3.6459 mL	7.2918 mL
50 mM	0.1458 mL	0.7292 mL	1.4584 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Wang J, Niu Y, Zhang C, Chen Y. A micro-plate colorimetric assay for rapid determination of trace zinc in animal feed, pet food and drinking water by ion masking and statistical partitioning correction. *Food Chem.* 2018 Apr 15; 245:337-345. doi: 10.1016/j.foodchem.2017.10.054. Epub 2017 Oct 14. PubMed PMID: 29287379.

McLellan R, Palacios MA, Sanz S, Brechin EK, Dalgarno SJ. Importance of Steric Influences in the Construction of Multicomponent Hybrid Polymetallic Clusters. *Inorg Chem.* 2017 Aug 21;56(16):10044-10053. doi: 10.1021/acs.inorgchem.7b01566. Epub 2017 Aug 3. PubMed PMID: 28771002.

Wei Z, Liu YQ, Wang SZ, Yao L, Nie HF, Wang YA, Liu XY, Zheng ZB, Li S. Conjugates of salicylaldoximes and peripheral site ligands: Novel efficient nonquaternary reactivators for nerve agent-inhibited acetylcholinesterase. *Bioorg Med Chem.* 2017 Aug 15;25(16):4497-4505. doi: 10.1016/j.bmc.2017.06.041. Epub 2017 Jun 27. PubMed PMID: 28684009.

Mills S E , Levine L . The Inhibition of Immune Haemolysis by Salicylaldoxime[J]. *Immunology*, 1959, 2(4):368-83.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481