

Forskolin

Chemical Properties

CAS No. : 66575-29-9

Formula: C₂₂H₃₄O₇


Molecular Weight: 410.50

Storage:

Keep away from direct sunlight, Keep away from moisture, Store at low temperature

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|---------------|---|
| Description | Forskolin (Coleonol) is a natural product, an adenylate cyclase activator (EC ₅₀ =0.5 μM). Forskolin increases cAMP levels, activates PXR and FXR, and induces autophagy. Forskolin produces positive inotropic effects in the heart, and has platelet anticoagulant and antihypertensive effects. |
| Targets(IC50) | FXR, AChR, Autophagy, PKC, Adenylate cyclase |
| In vitro | <p>METHODS: Rat adrenal medullary chromosomal tumor cells PC12 were treated with Forskolin (0.01-10 μM) for 3-48 h. The growth inhibition of the cells was detected by MTT.</p> <p>RESULTS: The cell viability decreased rapidly after treatment with 10 μM Forskolin (88.4% after 6 h and 60.5% after 48 h). [1]</p> <p>METHODS: Human myeloma cells U266, H929, INA-6, RPMI 8226 and OPM-2 were treated with Forskolin (1-100 μM) for 72 h. Cell death was detected by Flow Cytometry.</p> <p>RESULTS: Forskolin dose-dependently induced cell death in human myeloma cells, with U266, OPM-2 and INA-6 being more sensitive than H929 and RPMI 8226 cells. [2]</p> <p>METHODS: Human IL-2-dependent leukemia cells Kit 225 and human leukemia cells MT-2 were treated with Forskolin (1-100 μM) for 20 min, and cAMP concentration was measured by ELISA.</p> <p>RESULTS: Forskolin induced an up-regulation of cAMP levels, which reached a maximum</p> |
| In vivo | <p>METHODS: To assay anti-tumor activity in vivo, Forskolin (4-5 mg/kg in PBS/DMSO solution (15:1)) was injected intraperitoneally into BALB/c nude mice bearing murine multiple myeloma tumor MOPC315, and was administered on the 2nd/4th/6th day after tumor cell injection.</p> <p>RESULTS: All mice eventually developed tumors, but Forskolin significantly delayed tumor growth in vivo. Compounds that increase cAMP may have therapeutic potential in the treatment of multiple myeloma. [4]</p> <p>METHODS: To investigate the effect of Forskolin on retinal inflammation under diabetic conditions, Forskolin (50 mg/kg) was administered by gavage once weekly for twelve weeks to C57BL/6 mice in the STZ-induced diabetes model.</p> <p>RESULTS: Retinal glucose concentrations were increased in both diabetic control and Forskolin-treated groups compared to normal controls, but due to down-regulation of glucose transporter protein 1 expression, the Forskolin-treated group was only about</p> |

| | |
|-----------------|--|
| In vivo | 68.06% of the diabetic control group. ICAM-1 and TNF- α expression was up-regulated in the Forskolin-treated and diabetic control groups compared to the normal control group, but the expression levels of these two inflammatory factors in the Forskolin-treated group were 68.75% and 75.37% of those in the diabetic control group, respectively. [5] |
| Kinase Assay | For Jak3 kinase assays, Fsk-treated MT-2 cells were lysed, clarified, and immunoprecipitated using Jak3 antibody as described above. Kinase reactions were carried out as described previously at 30 °C for 20 min. For PKA kinase assays, untreated MT-2 cells were lysed, and Jak3 was immunoprecipitated and bound to PAS beads as described previously. Immunoprecipitated Jak3 was washed with kinase buffer (50 mM HEPES-NaOH (pH 7.4), 10 mM MgCl ₂ , 0.5 mM EGTA, 0.5 mM DTT, 20 μ g/ml aprotinin, 10 μ g/ml leupeptin, 1 μ g/ml pepstatin A) and incubated with 200 μ M ATP and purified protein kinase A catalytic subunit (PKAc) as indicated in the figure legends. Kinase reactions were carried out at 32 °C for 30 min followed by vigorous washing of the beads with cold kinase wash buffer as described previously. For [γ - ³² P]ATP radiolabeled kinase assays using recombinant Jak3, Hek293 cells were transfected with wild type (WT) Jak3 or kinase-dead Jak3 K855A using Lipofectamine 2000 according to the manufacturer's instructions. Cells were lysed and immunoprecipitated with Jak3 antibody. Jak3-bound PAS beads were washed three times in cold lysis buffer followed by kinase buffer. Kinase reactions were initiated by adding 10 μ Ci [γ - ³² P]ATP, 10 μ M unlabeled ATP, and 1 μ g of purified PKAc to Jak3-bound PAS bead reaction mixtures. Kinase reactions were performed at 32 °C for 30 min. Jak3-bound PAS beads were washed three times in radioimmunoassay buffer (10 mM Tris-HCl, pH 7.4, 75 mM NaCl, 20 mM EDTA, 10 mM EGTA, 20 mM Na ₄ P ₂ O ₇ , 50 mM NaF, 20 mM 2-glycerolphosphate, 1 mM p-nitrophenyl phosphate, 0.1% Triton X-100) and one time in kinase wash buffer. The reactions were stopped by adding 2 \times SDS-PAGE sample buffer followed by SDS-PAGE. Coomassie stainable Jak3 bands were excised from the PVDF membrane and subjected to phosphoamino acid analysis [2]. |
| Cell Research | Kit 225 or MT-2 cells were treated with 1, 5, 10, 20, 50, or 100 μ M Forskolin for 20 min at 37 °C. Cells were lysed and clarified by centrifugation, and the concentration of cAMP was detected by direct cAMP ELISA. Optical density was measured at 405 nm, and the concentration of intracellular cAMP was calculated using a weighted four parameter logistic curve according to the manufactures instructions [2]. |
| Animal Research | Forskolin was dissolved in dimethyl sulfoxide (DMSO) and injected intraperitoneally into neonatal mice at postnatal days 4 (P4) and 5 (P5). Mice injected with DMSO served as the controls. The treated mice were euthanized at P6, and their retinas were isolated for whole-mount immunohistochemistry (IHC). We first tested the effect of different concentrations of forskolin on the survival rate and retinal vasculature and determined the optimal concentration, 1.0 μ g/50 μ L (0.3 mg/kg) at P4 and 1.5 μ g/50 μ L (0.5 mg/kg) at P5, used to compare the retinal vascular phenotypes between WT mice and Mrp4-deficient mice [4]. . After acclimatization for 2 weeks, animals were randomly divided into four groups of eight rats each and treated for six consecutive weeks as follows: The first group was treated with CCl ₄ (50% CCl ₄ /corn oil; 0.5 mL \cdot kg ⁻¹ , i.p.) twice a week to induce liver fibrosis. The second group was given forskolin only at a dose of 10 mg \cdot kg ⁻¹ , i.p., dissolved in a DMSO/saline solution (1:49) five times a week. The third group was given both CCl ₄ and forskolin. The dose of forskolin used here was based on the results of our preliminary study. The fourth group served as the normal control, receiving vehicles only. At 24 h after the last injection, blood samples were collected from the retro-orbital plexus after light anesthesia with sodium pentobarbital (50 mg \cdot kg ⁻¹ , i.p.). Serum was separated by centrifugation at 3000 \times g for 10 min and was used for the assessment of liver functions. Rats were killed by cervical dislocation, and livers were removed and weighed. A portion of liver tissue was washed and homogenized to obtain |

A DRUG SCREENING EXPERT

| | |
|-----------------|---|
| Animal Research | a 20% (w-v?) homogenate, which was used for assessment of oxidative stress, inflammatory and fibrogenic markers. Another portion was placed in formalin for immunohistochemical and histopathological analyses. The remainder was stored at ?80° C, together with the 20% homogenate, until needed [5]. |
|-----------------|---|

Solubility Information

| | |
|---------------------|--|
| Solubility | H2O: Insoluble, DMSO: 245.00 mg/mL (596.83 mM),Sonication is recommended. Ethanol: 15.00 mg/mL (36.54 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble) |
| In vivo Formulation | 10% DMSO+40% PEG300+5% Tween 80+45% Saline: 3.00 mg/mL (7.31 mM),Solution. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i> |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 2.4361 mL | 12.1803 mL | 24.3605 mL |
| 5 mM | 0.4872 mL | 2.4361 mL | 4.8721 mL |
| 10 mM | 0.2436 mL | 1.218 mL | 2.4361 mL |
| 50 mM | 0.0487 mL | 0.2436 mL | 0.4872 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

- Park KH, et al. Modulation of PC12 cell viability by forskolin-induced cyclic AMP levels through ERK and JNK pathways: an implication for L-DOPA-induced cytotoxicity in nigrostriatal dopamine neurons. *Toxicol Sci.* 2012 Jul; 128(1):247-57.
- Li J, Bai Y, Liu Y, et al. Transcriptome-based chemical screens identify CDK8 as a common barrier in multiple cell reprogramming systems. *Cell Reports.* 2023, 42(6).
- Zhang X, Ke P X, Yuan X, et al. Forskolin Protected against Streptozotocin-Induced Diabetic Cardiomyopathy via Inhibition of Oxidative Stress and Cardiac Fibrosis in Mice. *BioMed Research International.* 2021, 2021
- Jin C, Zhao S, Xie H. Forskolin enhanced the osteogenic differentiation of human dental pulp stem cells in vitro and in vivo. *Journal of Dental Sciences.* 2022
- Wang H, Liang N, Huang D, et al. The effects of high-density lipoprotein and oxidized high-density lipoprotein on forskolin-induced syncytialization of BeWo cells. *Placenta.* 2021, 103: 199-205
- Miao Y, Zhang Y, Qiao S, et al. Oral administration of curcumin ameliorates pulmonary fibrosis in mice through 15d-PGJ2-mediated induction of hepatocyte growth factor in the colon. *Acta Pharmacologica Sinica.* 2020: 1-14
- Lin J Y, Cheng J, Du Y Q, et al. In vitro expansion of pancreatic islet clusters facilitated by hormones and chemicals. *Cell Discovery.* 2020, 6(1): 1-12
- Hu Z, Lu Y, Cao J, et al. N-acetyltransferase NAT10 controls cell fates via connecting mRNA cytidine acetylation to chromatin signaling. *Science Advances.* 2024, 10(2): eadh9871.
- Zhao X, Cai X, Zhu H, et al. 27-Hydroxycholesterol inhibits trophoblast fusion during placenta development by activating PI3K/AKT/mTOR signaling pathway. *Archives of Toxicology.* 2024: 1-15.
- Mendonça L S, Henriques D, Fernandes V, et al. Graft-derived neurons and bystander effects are maintained for six months after human iPSC-derived NESC transplantation in mice's cerebella. *Scientific Reports.* 2024, 14(1): 3236.
- Dang Q, Zhu Y, Zhang Y, et al. Nuclear Binding Protein 2/Nesfatin-1 Affects Trophoblast Cell Fusion during Placental Development via the EGFR-PLCG1-CAMK4 Pathway. *International Journal of Molecular Sciences.* 2024, 25 (3): 1925.
- Zhang Q, Shen X, Yuan X, et al. Lipopolysaccharide binding protein resists hepatic oxidative stress by regulating lipid droplet homeostasis. *Nature Communications.* 2024, 15(1): 3213.
- Follin-Arbelet V, et al. The natural compound forskolin synergizes with dexamethasone to induce cell death in myeloma cells via BIM. *Sci Rep.* 2015 Aug 26;5:13001.
- Han L, Song B, Zhang P, et al. PC3T: a signature-driven predictor of chemical compounds for cellular transition. *Communications Biology.* 2023, 6(1): 989.
- Rodriguez G, et al. Forskolin-inducible cAMP pathway negatively regulates T-cell proliferation by uncoupling the interleukin-2 receptor complex. *J Biol Chem.* 2013 Mar 8;288(10):7137-46.
- Chen S, Zhou X, Li W, et al. Development of a novel peptide targeting GPR81 to suppress adipocyte-mediated tumor progression. *Biochemical Pharmacology.* 2023: 115800.
- Ma X, Lu Y, Zhou Z, Human expandable pancreatic progenitor-derived β cells ameliorate diabetes. *Science Advances.* 2022, 8(8): eabk1826.
- Follin-Arbelet V, et al. Cyclic AMP induces apoptosis in multiple myeloma cells and inhibits tumor development in a mouse myeloma model. *BMC Cancer.* 2011 Jul 18;11:301.
- You ZP, et al. Forskolin attenuates retinal inflammation in diabetic mice. *Mol Med Rep.* 2018 Feb;17(2):2321-2326.
- Fu J, Jiang L, Yu B, et al. Generation of a Human iPSC Line CIBi010-A with a Reporter for ASGR1 Using CRISPR/Cas9. *Stem Cell Research.* 2022: 102800
- Lu J, Dou F, Yu Z. The potassium channel KCa3. 1 represents a valid pharmacological target for microgliosis-induced neuronal impairment in a mouse model of Parkinson's disease[J]. *Journal of Neuroinflammation.* 2019, 16 (1): 1-14.
- Zhao C, Sun C, Yuan J, et al. *Heridium caput-medusae* (Bull.: Fr.) Pers. fermentation concentrate polysaccharides improves intestinal bacteria by activating chloride channels and mucus secretion. *Journal of Ethnopharmacology.* 2022: 115721.
- Xiaoli F, Yaqing Z, Ruhui L, et al. Graphene oxide disrupted mitochondrial homeostasis through inducing intracellular redox deviation and autophagy-lysosomal network dysfunction in SH-SY5Y cells. *Journal of Hazardous Materials.* 2021: 126158.
- Lin J Y, Cheng J, Du Y Q, et al. In vitro expansion of pancreatic islet clusters facilitated by hormones and chemicals

]]. Cell Discovery. 2020, 6(1): 1-12.

Wu L, Dong A, Dong L, et al. PARIS, an optogenetic method for functionally mapping gap junctions]]. eLife. 2019 Jan 14;8. pii: e43366.

Wu L, Dong A, Dong L, et al. PARIS, an optogenetic method for functionally mapping gap junctions. ELife. 2019 Jan 14;8. pii: e43366

Wang W, Ren S, Lu Y, et al. Inhibition of Syk promotes chemical reprogramming of fibroblasts via metabolic rewiring and H2S production. The EMBO Journal. 2021 Jun 1;40(11):e106771. doi: 10.15252/embj.2020106771.

Epub 2021 Apr 28.

Xiaoli F, Yaqing Z, Ruhui L, et al. Graphene oxide disrupted mitochondrial homeostasis through inducing intracellular redox deviation and autophagy-lysosomal network dysfunction in SH-SY5Y cells]]. Journal of Hazardous Materials. 2021: 126158.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481