


L162441

Chemical Properties

CAS No. : 154512-46-6
 Formula: C₃₄H₄₃N₅O₅S
 Molecular Weight: 633.8
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	L162441 is an antagonist of Angiotensin type 1 receptor.
Targets(IC50)	RAAS

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5778 mL	7.8889 mL	15.7778 mL
5 mM	0.3156 mL	1.5778 mL	3.1556 mL
10 mM	0.1578 mL	0.7889 mL	1.5778 mL
50 mM	0.0316 mL	0.1578 mL	0.3156 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481