


4-Chloro-1H-pyrazolo[3,4-d]pyrimidine

Chemical Properties

CAS No. :	5399-92-8
Formula:	C ₅ H ₃ ClN ₄
Molecular Weight:	154.56
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Heterocyclic compound-purine, intermediate and building block-electrophile, nucleoside base
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.470 mL	32.3499 mL	64.6998 mL
5 mM	1.294 mL	6.470 mL	12.940 mL
10 mM	0.647 mL	3.235 mL	6.470 mL
50 mM	0.1294 mL	0.647 mL	1.294 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481