


SCS

## Chemical Properties

CAS No. : 3232-36-8  
 Formula: C<sub>14</sub>H<sub>12</sub>N<sub>2</sub>O<sub>3</sub>  
 Molecular Weight: 256.26  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



## Biological Description

Description	SCS is a GABAA receptor antagonist.
Targets(IC50)	GABA Receptor

## Solubility Information

Solubility	DMSO: 50 mg/mL (195.11 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween-80+45% Saline: 1 mg/mL (3.9 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9023 mL	19.5114 mL	39.0229 mL
5 mM	0.7805 mL	3.9023 mL	7.8046 mL
10 mM	0.3902 mL	1.9511 mL	3.9023 mL
50 mM	0.078 mL	0.3902 mL	0.7805 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Thompson SA, et al. Salicylidene salicylhydrazide, a selective inhibitor of beta 1-containing GABAA receptors. Br J Pharmacol. 2004 May;142(1):97-106.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481