

Estril 3-sulfate

Chemical Properties


CAS No. :

Formula: C₁₈H₂₄O₆S

Molecular Weight: 368.44

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Estril 3-sulfate is a useful organic compound for research related to life sciences and the catalog number is T124799.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7141 mL	13.5707 mL	27.1415 mL
5 mM	0.5428 mL	2.7141 mL	5.4283 mL
10 mM	0.2714 mL	1.3571 mL	2.7141 mL
50 mM	0.0543 mL	0.2714 mL	0.5428 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481