


Alpibectir

Chemical Properties

CAS No. : 2285440-39-1
 Formula: C₁₂H₁₄F₆N₂O₂
 Molecular Weight: 332.24
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Alpibectir has antibacterial activity [1].
Targets(IC50)	Others,Antibacterial

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0099 mL	15.0494 mL	30.0987 mL
5 mM	0.602 mL	3.0099 mL	6.0197 mL
10 mM	0.301 mL	1.5049 mL	3.0099 mL
50 mM	0.0602 mL	0.301 mL	0.602 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481