


## (2-Chloro-4-pyridinyl)methanol

## Chemical Properties

CAS No. :	100704-10-7
Formula:	C <sub>6</sub> H <sub>6</sub> ClNO
Molecular Weight:	143.57
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(2-Chloro-4-pyridinyl)methanol, with CAS No. 100704-10-7, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. (2-Chloro-4-pyridinyl)methanol provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.9652 mL	34.8262 mL	69.6524 mL
5 mM	1.393 mL	6.9652 mL	13.9305 mL
10 mM	0.6965 mL	3.4826 mL	6.9652 mL
50 mM	0.1393 mL	0.6965 mL	1.393 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481