

cIAP1 Ligand-Linker Conjugates 2 Hydrochloride

Chemical Properties

CAS No. :


Formula: C₃₇H₄₉ClN₄O₇

Molecular Weight: 697.26

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	cIAP1 Ligand-Linker Conjugates 2 Hydrochloride is a chemical compound that combines an IAP ligand, targeting the E3 ubiquitin ligase, with a PROTAC linker and is utilized in the creation of SNIPERs[1].
Targets(IC50)	Others,E3 Ligase Ligand-Linker Conjugates,PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4342 mL	7.1709 mL	14.3419 mL
5 mM	0.2868 mL	1.4342 mL	2.8684 mL
10 mM	0.1434 mL	0.7171 mL	1.4342 mL
50 mM	0.0287 mL	0.1434 mL	0.2868 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Itoh Y, et al. Design, synthesis and biological evaluation of nuclear receptor-degradation inducers. Bioorg Med Chem. 2011 Nov 15;19(22):6768-78.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481