

MS 39N

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	MS 39N is a negative control for MS 39. Binds EGFR without inducing degradation.
Targets(IC50)	Others

Reference

Cheng et al (2020) Discovery of potent and selective epidermal growth factor receptor (EGFR) bifunctional small-molecule degraders. J.Med.Chem. 63 1216 PMID: 31895569

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481