


Dihydropteroate synthase-IN-1

Chemical Properties

CAS No. : 2418026-70-5
 Formula: C₁₉H₂₃N₅O₄S₂
 Molecular Weight: 449.55
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Dihydropteroate synthase-IN-1 (compound 5g) is a potent inhibitor of dihydropteroate synthase (DHPS) and cytochrome P450, and can be used as a diagnostic radiographic material.
Targets(IC50)	Others,Antibacterial

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2244 mL	11.1222 mL	22.2445 mL
5 mM	0.4449 mL	2.2244 mL	4.4489 mL
10 mM	0.2224 mL	1.1122 mL	2.2244 mL
50 mM	0.0445 mL	0.2224 mL	0.4449 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481