

RO5073012

## Chemical Properties


CAS No. : 1020814-70-3

Formula: C<sub>13</sub>H<sub>16</sub>ClN<sub>3</sub>

Molecular Weight: 249.74

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	RO5073012 is a novel potent and selective TAAR1 partial agonist.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.0042 mL	20.0208 mL	40.0416 mL
5 mM	0.8008 mL	4.0042 mL	8.0083 mL
10 mM	0.4004 mL	2.0021 mL	4.0042 mL
50 mM	0.0801 mL	0.4004 mL	0.8008 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481