

Propargyl-C2-NHS ester

Chemical Properties

CAS No. :	906564-59-8
Formula:	C ₁₀ H ₁₁ NO ₄
Molecular Weight:	209.199
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Propargyl-C2-NHS ester is a noncleavable linker used in antibody-drug conjugation (ADC).
Targets(IC50)	ADC Linker,PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.7801 mL	23.9006 mL	47.8011 mL
5 mM	0.956 mL	4.7801 mL	9.5602 mL
10 mM	0.478 mL	2.3901 mL	4.7801 mL
50 mM	0.0956 mL	0.478 mL	0.956 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481