


## Yunaconitoline

## Chemical Properties

CAS No. :	259099-25-7
Formula:	C <sub>34</sub> H <sub>43</sub> NO <sub>10</sub>
Molecular Weight:	625.71
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Yunaconitoline is a natural product for research related to life sciences. The catalog number is TN5273 and the CAS number is 259099-25-7.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5982 mL	7.9909 mL	15.9818 mL
5 mM	0.3196 mL	1.5982 mL	3.1964 mL
10 mM	0.1598 mL	0.7991 mL	1.5982 mL
50 mM	0.032 mL	0.1598 mL	0.3196 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481