


AMG0074

Chemical Properties

CAS No. : 871507-11-8
 Formula: C₂₄H₁₈N₄O
 Molecular Weight: 378.43
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	AMG0074 is a PDE10 antagonist.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6425 mL	13.2125 mL	26.425 mL
5 mM	0.5285 mL	2.6425 mL	5.285 mL
10 mM	0.2642 mL	1.3212 mL	2.6425 mL
50 mM	0.0528 mL	0.2642 mL	0.5285 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481