


Picilorex

Chemical Properties

CAS No. :	62510-56-9
Formula:	C ₁₄ H ₁₈ ClN
Molecular Weight:	235.75
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Picilorex is a monoamine reuptake inhibitor, a stimulant as well as a derivate of Pyrrolidine.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.2418 mL	21.2089 mL	42.4178 mL
5 mM	0.8484 mL	4.2418 mL	8.4836 mL
10 mM	0.4242 mL	2.1209 mL	4.2418 mL
50 mM	0.0848 mL	0.4242 mL	0.8484 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481