


3'-Azido-3'-deoxy-beta-L-adenosine

Chemical Properties

CAS No. :	216976-76-0
Formula:	C10H12N8O3
Molecular Weight:	292.25
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3'-Azido-3'-deoxy-beta-L-adenosine is a Nucleoside Derivative - L-nucleoside, Azido-nucleoside.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4217 mL	17.1086 mL	34.2173 mL
5 mM	0.6843 mL	3.4217 mL	6.8435 mL
10 mM	0.3422 mL	1.7109 mL	3.4217 mL
50 mM	0.0684 mL	0.3422 mL	0.6843 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481