

## PROTAC IRAK4 degrader-5

## Chemical Properties

CAS No. : 2360530-61-4


Formula: C41H40F3N11O9

Molecular Weight: 887.834

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	PROTAC IRAK4 degrader-5 is a Cereblon-based compound designed to selectively degrade IRAK4.
Targets(IC50)	Others,IRAK,PROTACs

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1263 mL	5.6317 mL	11.2634 mL
5 mM	0.2253 mL	1.1263 mL	2.2527 mL
10 mM	0.1126 mL	0.5632 mL	1.1263 mL
50 mM	0.0225 mL	0.1126 mL	0.2253 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Mainolfi N, et, al. Irak degraders and uses thereof. US20190192668A1.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481