

NPD8716

## Chemical Properties


CAS No. : 956944-02-8

Formula: C<sub>23</sub>H<sub>31</sub>N<sub>1</sub>O<sub>6</sub>S

Molecular Weight: 449.56

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	NPD8716 is a hepatitis B/D virus entry inhibitor, targeting sodium taurocholate cotransporting polypeptide (NTCP).
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2244 mL	11.122 mL	22.244 mL
5 mM	0.4449 mL	2.2244 mL	4.4488 mL
10 mM	0.2224 mL	1.1122 mL	2.2244 mL
50 mM	0.0445 mL	0.2224 mL	0.4449 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481