


2,3-dihydroxy-3-(4-hydroxyphenyl)propanoic acid

Chemical Properties

CAS No. :	100201-57-8
Formula:	C ₉ H ₁₀ O ₅
Molecular Weight:	198.17
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2,3-dihydroxy-3-(4-hydroxyphenyl)propanoic acid is a natural product for research related to life sciences. The catalog number is TN2702 and the CAS number is 100201-57-8.
Targets(IC50)	Others,Endogenous Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.0462 mL	25.2309 mL	50.4617 mL
5 mM	1.0092 mL	5.0462 mL	10.0923 mL
10 mM	0.5046 mL	2.5231 mL	5.0462 mL
50 mM	0.1009 mL	0.5046 mL	1.0092 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481