


3'-O-(t-Butyldiphenylsilyl)thymidine

Chemical Properties

CAS No. :	83467-48-5
Formula:	C ₂₆ H ₃₂ N ₂ O ₅ Si
Molecular Weight:	480.63
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3'-O-(t-Butyldiphenylsilyl)thymidine is a Nucleoside Derivative - Protected nucleoside w/NH ₂ /OH open.
Targets(IC ₅₀)	Nucleoside Antimetabolite/Analog,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0806 mL	10.403 mL	20.806 mL
5 mM	0.4161 mL	2.0806 mL	4.1612 mL
10 mM	0.2081 mL	1.0403 mL	2.0806 mL
50 mM	0.0416 mL	0.2081 mL	0.4161 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481