


Senkyunolide I

Chemical Properties

CAS No. :	94596-28-8
Formula:	C ₁₂ H ₁₆ O ₄
Molecular Weight:	224.25
Storage:	Pure form: -20°C for 3 years In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	1. Senkyunolide I has neuroprotective effect, associated with its anti-oxidation and anti-apoptosis properties.
Targets(IC50)	Caspase

Solubility Information

Solubility	DMSO: 55 mg/mL (245.26 mM),Sonication is recommended. H2O: 40 mg/mL (178.37 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (8.92 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.4593 mL	22.2965 mL	44.5931 mL
5 mM	0.8919 mL	4.4593 mL	8.9186 mL
10 mM	0.4459 mL	2.2297 mL	4.4593 mL
50 mM	0.0892 mL	0.4459 mL	0.8919 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Hu Y , Duan M , Liang S , et al. Senkyunolide I protects rat brain against focal cerebral ischemia-reperfusion injury by up-regulating p-Erk1/2, Nrf2/HO-1 and inhibiting caspase 3[J]. Brain Research, 2015, 1605:39-48.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481