


## Neostigmine methyl sulfate

## Chemical Properties

CAS No. :	51-60-5
Formula:	C <sub>13</sub> H <sub>22</sub> N <sub>2</sub> O <sub>6</sub> S
Molecular Weight:	334.39
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Neostigmine methyl sulfate (Syntostigmin) is a reversible inhibitor of acetylcholinesterase, can not cross the blood-brain barrier, used in the treatment of myasthenia gravis.
Targets(IC50)	AChR, Cholinesterase (ChE)

## Solubility Information

Solubility	DMSO: 55 mg/mL (164.48 mM), Sonication is recommended. H <sub>2</sub> O: 200.4 mM, Sonication is recommended. ( $< 1$ mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (5.98 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9905 mL	14.9526 mL	29.9052 mL
5 mM	0.5981 mL	2.9905 mL	5.981 mL
10 mM	0.2991 mL	1.4953 mL	2.9905 mL
50 mM	0.0598 mL	0.2991 mL	0.5981 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Huang TH, et al. Yao Xue Xue Bao. 1982 Jul;17(7):534-9.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481