


LysoFP-NH2

Chemical Properties

CAS No. :	69408-85-1
Formula:	C18H19N3O3
Molecular Weight:	325.368
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	LysoFP-NH2 is the fluorescent form of the lysosomal turn-on fluorescent probe for carbon monoxide (CO) lysoFP-NO2. LysoFP-NH2 is formed when lysoFP-NO2 reacts with CO. It displays excitation/emission maxima of 440/528 nm, respectively.
Targets(IC50)	Others

Solubility Information

Solubility	Ethanol: 1 mg/mL (3.07 mM), Sonication is recommended. DMSO:PBS (pH 7.2) (1:30): 0.03 mg/mL (0.09 mM), Sonication is recommended. DMSO: 20 mg/mL (61.47 mM), Sonication is recommended. DMF: 20 mg/mL (61.47 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0734 mL	15.3671 mL	30.7342 mL
5 mM	0.6147 mL	3.0734 mL	6.1468 mL
10 mM	0.3073 mL	1.5367 mL	3.0734 mL
50 mM	0.0615 mL	0.3073 mL	0.6147 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481