


YD2

Chemical Properties

CAS No. : 1201915-91-4
 Formula: C₈₁H₈₉N₅O₂Zn
 Molecular Weight: 1230.02
 Storage: Keep away from direct sunlight
 Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	YD2 is a porphyrin sensitizer. YD2 is a useful material for making solar cell. Cosensitization of YD2-o-C8 with another organic dye enhances the performance of the device, leading to a measured power conversion efficiency of 12.3% under simulated air mass 1.5 global sunlight.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.813 mL	4.065 mL	8.1299 mL
5 mM	0.1626 mL	0.813 mL	1.626 mL
10 mM	0.0813 mL	0.4065 mL	0.813 mL
50 mM	0.0163 mL	0.0813 mL	0.1626 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Chiba S, Wakabayashi T, Tachi Y, Takaki N, Terashima A, Okumura S, Yoshida T. Method to Reduce Long-lived Fission Products by Nuclear Transmutations with Fast Spectrum Reactors. *Sci Rep.* 2017 Oct 24;7(1):13961. doi: 10.1038/s41598-017-14319-7. Erratum in: *Sci Rep.* 2017 Dec 5;7(1):17270. PubMed PMID: 29066843; PubMed Central PMCID: PMC5654822.

Balanay MP, Kim DH. Strategic Design of Bacteriochlorins as Possible Dyes for Photovoltaic Applications. *J Phys Chem A.* 2017 Sep 7;121(35):6660-6669. doi: 10.1021/acs.jpca.7b07957. Epub 2017 Aug 28. PubMed PMID: 28800232.

Higashino T, Kurumisawa Y, Cai N, Fujimori Y, Tsuji Y, Nimura S, Packwood DM, Park J, Imahori H. A Hydroxamic Acid Anchoring Group for Durable Dye-Sensitized Solar Cells Incorporating a Cobalt Redox Shuttle. *ChemSusChem.* 2017 Sep 11;10(17):3347-3351. doi: 10.1002/cssc.201701157. Epub 2017 Aug 16. PubMed PMID: 28722326.

Zhang CR, Li XY, Shen YL, Wu YZ, Liu ZJ, Chen HS. Molecular Docking toward Panchromatic Dye Sensitizers for Solar Cells Based upon Tetraazulenylporphyrin and Tetraanthracenylporphyrin. *J Phys Chem A.* 2017 Apr 6;121(13):2655-2664. doi: 10.1021/acs.jpca.6b12979. Epub 2017 Mar 28. PubMed PMID: 28319383.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481