

## Compound T127974

## Chemical Properties

CAS No. :

Formula: C<sub>26</sub>H<sub>42</sub>O<sub>10</sub>

Molecular Weight: 514.61

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	4b-hydroxy-1,4a,7-trimethyl-10-oxo-7-(2-([3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy)ethyl)-tetradecahydrophenanthrene-1-carboxylic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9432 mL	9.7161 mL	19.4322 mL
5 mM	0.3886 mL	1.9432 mL	3.8864 mL
10 mM	0.1943 mL	0.9716 mL	1.9432 mL
50 mM	0.0389 mL	0.1943 mL	0.3886 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481