


## C.I. Vat Black 27

## Chemical Properties

CAS No. :	2379-81-9
Formula:	C42H23N3O6
Molecular Weight:	665.65
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	C.I. Vat Black 27 is a black dye.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5023 mL	7.5115 mL	15.0229 mL
5 mM	0.3005 mL	1.5023 mL	3.0046 mL
10 mM	0.1502 mL	0.7511 mL	1.5023 mL
50 mM	0.030 mL	0.1502 mL	0.3005 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481