


Diosmetinidin chloride

Chemical Properties

CAS No. :	64670-94-6
Formula:	C ₁₆ H ₁₃ ClO ₅
Molecular Weight:	320.72
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Diosmetinidin chloride is a natural product for research related to life sciences. The catalog number is TN6616 and the CAS number is 64670-94-6.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.118 mL	15.5899 mL	31.1798 mL
5 mM	0.6236 mL	3.118 mL	6.236 mL
10 mM	0.3118 mL	1.559 mL	3.118 mL
50 mM	0.0624 mL	0.3118 mL	0.6236 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481