


Pentopril

Chemical Properties

CAS No. : 82924-03-6
 Formula: C₁₈H₂₃NO₅
 Molecular Weight: 333.38
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Pentopril is an Enzyme Inhibitor (Angiotensin-Converting)
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9996 mL	14.9979 mL	29.9958 mL
5 mM	0.5999 mL	2.9996 mL	5.9992 mL
10 mM	0.300 mL	1.4998 mL	2.9996 mL
50 mM	0.060 mL	0.300 mL	0.5999 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481