


Dehydrogeijerin

Chemical Properties

CAS No. :	16850-91-2
Formula:	C ₁₅ H ₁₄ O ₄
Molecular Weight:	258.27
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Dehydrogeijerin significantly inhibits the inflammatory activity of activated macrophages, suggesting that it could be a potential candidate for the treatment of inflammatory disease.
Targets(IC50)	NOS,MAPK,Cholinesterase (ChE),COX

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8719 mL	19.3596 mL	38.7192 mL
5 mM	0.7744 mL	3.8719 mL	7.7438 mL
10 mM	0.3872 mL	1.936 mL	3.8719 mL
50 mM	0.0774 mL	0.3872 mL	0.7744 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481