

## Elgonica dimer A

## Chemical Properties


CAS No. :

Formula: C<sub>36</sub>H<sub>30</sub>O<sub>14</sub>

Molecular Weight: 686.622

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Elgonica dimer A is a useful organic compound for research related to life sciences and the catalog number is T126051.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4564 mL	7.282 mL	14.5641 mL
5 mM	0.2913 mL	1.4564 mL	2.9128 mL
10 mM	0.1456 mL	0.7282 mL	1.4564 mL
50 mM	0.0291 mL	0.1456 mL	0.2913 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481