


Taraxinic acid

Chemical Properties

CAS No. :	75911-33-0
Formula:	C ₁₅ H ₁₈ O ₄
Molecular Weight:	262.305
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Taraxinic acid is a natural product for research related to life sciences. The catalog number is TN6236 and the CAS number is 75911-33-0.
Targets(IC50)	Tyrosinase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8123 mL	19.0614 mL	38.1228 mL
5 mM	0.7625 mL	3.8123 mL	7.6246 mL
10 mM	0.3812 mL	1.9061 mL	3.8123 mL
50 mM	0.0762 mL	0.3812 mL	0.7625 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481