

Compound T127878

Chemical Properties


CAS No. :

Formula: C₃₂H₅₂O₅

Molecular Weight: 516.76

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

(1R,7R,8R,9aS,10R,11R,11aR)-7,8,10-trihydroxy-6,6,9a,11a-tetramethyl-1-(6-methyl-5-methylideneheptan-2-yl)-1H,2H,3H,3aH,4H,5H,5aH,6H,7H,8H,9H,9aH,10H,11H,11aH-cyclopenta[a]phenanthren-11-yl acetate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9351 mL	9.6757 mL	19.3513 mL
5 mM	0.387 mL	1.9351 mL	3.8703 mL
10 mM	0.1935 mL	0.9676 mL	1.9351 mL
50 mM	0.0387 mL	0.1935 mL	0.387 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481