


(S)-Pomalidomide

Chemical Properties

CAS No. :	202271-89-4
Formula:	C ₁₃ H ₁₁ N ₃ O ₄
Molecular Weight:	273.24
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(S)-Pomalidomide is an inhibitor of angiogenesis and growth of B-cell neoplasias, which induces complete regressions of Burkitt's lymphoma cell tumors.
Targets(IC50)	Others,Endogenous Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6598 mL	18.2989 mL	36.5979 mL
5 mM	0.732 mL	3.6598 mL	7.3196 mL
10 mM	0.366 mL	1.8299 mL	3.6598 mL
50 mM	0.0732 mL	0.366 mL	0.732 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481