


## 1-Phenylpiperazine

## Chemical Properties

CAS No. :	92-54-6
Formula:	C10H14N2
Molecular Weight:	162.23
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	1-Phenylpiperazine, with CAS No. 92-54-6, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 1-Phenylpiperazine provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.1641 mL	30.8204 mL	61.6409 mL
5 mM	1.2328 mL	6.1641 mL	12.3282 mL
10 mM	0.6164 mL	3.082 mL	6.1641 mL
50 mM	0.1233 mL	0.6164 mL	1.2328 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481